



**Wellell Group**

# PRODUCT INFORMATION

**universal hoists and active hoists**  
slings and accessories

#workgreen

[www.slk-gmbh.de](http://www.slk-gmbh.de)



# content

## **4/5 universal hoists**

- 6** SLK Carry 185 Eco
- 7** SLK Carry 185 Classic / Classic e
- 8** SLK Carry 185 Pro / Pro e
- 9** SLK Carry 185 Pro L / Pro L e
- 10** SLK Carry Compact

## **12/13 hybrid hoists**

- 14** SLK Multy Universal / Universal e
- 15** SLK Multy Conversion Kit Active
- 16** SLK Multy Active / Active e
- 17** SLK Multy Conversion Kit Universal

## **18 active hoists**

- 20** SLK Eazy-up

## **21 technical specifications universal hoists & active hoists**

- 22/23** Technical Specifications Universal Hoists & Active Hoists

## **24/25 universal slings / active slings**

- 26** SLK Universal Sling
- 27** SLK Universal Sling with Headrest
- 28** SLK Universal Bathroom Sling
- 29** SLK Universal Bathroom Sling with Headrest
- 30** SLK Quick Transfer and Toilet Sling
- 31** SLK Toilet Sling
- 32** SLK Hammock Universal Sling
- 33** SLK Hammock Bathroom Sling
- 34** SLK Hammock Seat Pan Sling
- 35** SLK Comfort Lying Sling
- 36** SLK Comfort Lying Bathroom Sling
- 37** SLK Back-Fix Sling
- 38** SLK Back-Fix Bathroom Sling
- 39** SLK Maxi Sling
- 40** SLK Active Sling
- 41** SLK Special Active Sling
- 42** SLK Special Active Sling with Seat Swing

## **43** size chart / overview

**44** Size Chart Slings

**45** Overview Slings

## **46** hoist accessories / slings accessories

**47** SLK Hoist Scale

**47** SLK Calibrated Hoist Scale

**48** SLK Scales incl. Scale Adapter

**48** SLK Charging Station

**48** SLK Charging Station with Charger

**48** SLK Wall-mounted Charging Station

**49** SLK 2-Point Spreader Bar

**49** SLK 4-Point Spreader Bar X

**49** SLK 4-Point Spreader Bar H

**50** SLK 8-Point Spreader Bar

**51** SLK Impact Protector for Lifting Arm

**51** SLK Velcro Headrest

**51** SLK Adjustable Headrest Reinforced

**52** SLK Extension Loops

**53** SLK Weight Test Trolley up to 220 kg

**54** SLK Test Trolley up to 440 kg

# #workgreen

# For these reasons you should use Universal Hoists

## How many people are needed according to load handling regulation to reposition a patient?

According to national law, 4 nurses are required. Alternatively, a patient hoist can be used which is operated by only one nurse.

Universal Hoists have been used for decades in everyday nursing care to support care. They can be used for home care, in care facilities and hospitals. Universal Hoists are used to pick up, transport and change the position of patients with mobility impairments due to illness or disability. With a suitable Sling, it is possible to move a patient safely and gently.

## Advantages of a Universal Hoist

### Advantages for the patient

- safe and gentle transfer
- gentle pick up and settle of the patient
- lifting from floor level possible
- electric spreading (optional) adjustable

### Advantages for the caregiver

- simple operation during transfer
- minimization of physical stress on the caregiver during transfer

universal hoists  
**designed for use in facilities,  
clinics and homecare**



# SLK Carry 185 Eco

economical Universal Hoist for home care



DATA	MANUAL SPREADING
length	118 cm (46.4")
total height	137 cm (54")
width	61 cm (24")
leg height	11 cm (4.3")
lifting height	81 - 203 cm (32" - 80")
turning radius	132 cm (52")
base expansion	51 - 94 cm (20" - 37")
weight	37 kg (81 lbs)

further details on page: 22-23

### facts

- max. patient weight: 185 kg (407 lbs)
- protects the back of the caregiver
- gentle pick up and settle of the patient
- lifting from floor level possible
- manual spreading
- high manoeuvrability in a small area
- Hook mount with spring safety device
- 360-degree rotation of the Spreader Bar in highest possible position
- foldable without tools / transportable

### appropriate slings

- all SLK Slings (from page 25)

### standard

- SLK 2-Point Spreader Bar

### optional

- SLK 4-Point Spreader Bar X: 8060
- SLK 4-Point Spreader Bar H: 8064

Accessories from page 46

### material

- powder coated steel
- ball-bearing swivel-castors

### warranty

- According to our current service guidelines

medical aid number: (Germany) 22.40.01.0098 (manual)  
article-ID 8000 manual spreading

# SLK Carry 185 Classic / Classic e

ideal solution for home care



DATA	MANUAL SPREADING	ELECTRIC SPREADING
length	118 cm (46.4")	118 cm (46.4")
total height	137 cm (54")	137 cm (54")
width	61 cm (24")	61 cm (24")
leg height	11 cm (4.3")	11 cm (4.3")
lifting height	81 - 203 cm (31.8" - 80")	81 - 203 cm (32" - 80")
turning radius	132 cm (52")	132 cm (52")
base expansion	51 - 94 cm (20" - 37")	51 - 94 cm (20" - 37")
weight	37 kg (81 lbs)	38 kg (84 lbs)

further details on page: 22-23

## facts

- max. patient weight: 185 kg (407 lbs)
- protects the back of the caregiver
- gentle pick up and settle of the patient
- lifting from floor level possible
- electric spreading (optional) adjustable
- high manoeuvrability in a small area
- Hook mount with spring safety device
- 360-degree rotation of the Spreader Bar in highest possible position
- comfort-control with display and Readout function
- handling via control box possible
- soft-start / soft-stop
- high battery power
- foldable without tools / transportable
- standard mast and chassis are removable without tools

## appropriate slings

- all SLK Slings (from page 25)

## standard

- SLK 2-Point Spreader Bar

## optional

- SLK 4-Point Spreader Bar X: 8060
- SLK 4-Point Spreader Bar H: 8064

Accessories from page 46

## material

- powder coated steel
- ball-bearing swivel-casters

## warranty


- According to our current service guidelines

**ALSO AVAILABLE WITH e-CHASSIS...**

medical aid number: (Germany) 22.40.01.0099 (manual), 22.40.01.0110 (electric)  
article-ID 8002 manual spreading, 8003 electric spreading

# SLK Carry 185 Pro / Pro e

Universal Hoist designed for use in facilities and clinics



Chassis with comfort castors

DATA	MANUAL SPREADING	ELECTRIC SPREADING
length	115 cm (45.3")	115 cm (45.3")
total height	137 cm (54")	137 cm (54")
width	64 cm (26")	64 cm (26")
leg height	11 cm (4.3")	11 cm (4.3")
lifting height	81 - 203 cm (32" - 80")	81 - 203 cm (32" - 80")
turning radius	130 cm (51")	130 cm (51")
base expansion	51 - 94 cm (20" - 37")	51 - 94 cm (20" - 37")
weight	38 kg (84 lbs)	39 kg (86 lbs)

Further details on page: 22-23

## facts

- max. patient weight: 185 kg (407 lbs)
- protects the back of the caregiver
- gentle pick up and settle of the patient
- lifting from floor level possible
- electric spreading (optional) adjustable
- high manoeuvrability in a small area
- Hook mount with spring safety device
- 360-degree rotation of the Spreader Bar in highest possible position
- comfort-control with display and Readout function
- handling via control box possible
- soft-start / soft-stop
- high battery power
- foldable without tools / transportable
- standard mast and chassis are removable without tools

## appropriate slings

- all SLK Slings (from page 25)

## standard

- SLK 2-Point Spreader Bar

## optional

- SLK 4-Point Spreader Bar X: 8060
- SLK 4-Point Spreader Bar H: 8064

Accessories from page 46

## material

- powder coated steel
- ball-bearing, smooth-running comfort-castors with high-quality plastic housing

## warranty

- According to our current service guidelines


**ALSO AVAILABLE WITH e-CHASSIS...**

medical aid number: (Germany) 22.40.01.0100 (manual), 22.40.01.0108 (electric)  
article-ID 8004 manual spreading, 8005 electric spreading



# SLK Carry 185 Pro L / Pro L e

Universal Hoist with extra low chassis



Especially low chassis with comfort castors

DATA	MANUAL SPREADING	ELECTRIC SPREADING
length	115 cm (45.3")	115 cm (45.3")
total height	135 cm (53")	135 cm (53")
width	64 cm (26")	64 cm (26")
leg height	9 cm (3.5")	9 cm (3.5")
lifting height	79 - 201 cm (31" - 79")	79 - 201 cm (31" - 79")
turning radius	130 cm (51")	130 cm (51")
base expansion	51 - 94 cm (20" - 37")	51 - 94 cm (20" - 37")
weight	38 kg (84 lbs)	39 kg (86 lbs)

Further details on page: 22-23

## facts

- max. patient weight: 185 kg (407 lbs)
- protects the back of the caregiver
- gentle pick up and settle of the patient
- lifting from floor level possible
- electric spreading (optional) adjustable
- high manoeuvrability in a small area
- particularly low front underride height of: 88 mm
- Hook mount with spring safety device
- 360-degree rotation of the Spreader Bar in highest possible position
- comfort-control with display and Readout function
- handling via control box possible
- soft-start / soft-stop
- high battery power
- foldable without tools / transportable
- standard mast and chassis are removable without tools

## appropriate slings

- all SLK Slings (from page 25)

## standard

- SLK 2-Point Spreader Bar

## optional

- SLK 4-Point Spreader Bar X: 8060
- SLK 4-Point Spreader Bar H: 8064

Accessories from page 46

## material

- powder coated steel
- ball-bearing, smooth-running comfort-castors with high-quality plastic housing

## warranty

- According to our current service guidelines

**ALSO AVAILABLE WITH e-CHASSIS...**

medical aid number: (Germany) 22.40.01.0101 (manual), 22.40.01.0109 (electric)  
article-ID 8006 manual spreading, 8007 electric spreading

# SLK Carry Compact

Universal Hoist with short chassis for narrow living spaces



especially compact universal hoist

DATA	MANUAL SPREADING
length	103 cm (40.5")
total height	124 cm (50")
width	56 cm (22")
leg height	8,5 cm (3.4")
lifting height	81 - 175 cm (32" - 69")
turning radius	115 cm (45")
base expansion	45 - 78 cm (17.7" - 30.7")
weight	33 kg (73 lbs)

further details on page: 22-23

## facts

- max. patient weight: 160 kg (352 lbs)
- for children and small patients up to 1,70 m (67")
- protects the back of the caregiver
- gentle pick up and settle of the patient
- lifting from floor level possible
- adjustable chassis spreading
- easy to use in narrow spaces and doors
- high manoeuvrability in a small area
- particularly low front underride height of: 83 mm
- Hook mount with spring safety device
- 360-degree rotation of the Spreader Bar in highest possible position
- comfort-control with display and Readout function
- handling via control box possible
- soft-start / soft-stop
- high battery power / Li-Ion battery
- foldable without tools / transportable
- standard mast and chassis are removable without tools

## appropriate slings

- all SLK Slings (from page 25)

## standard

- SLK 2-Point Spreader Bar Compact H

## optional

- SLK 4-Point Spreader Bar H: 8064

Accessories from page 46

## material

- powder coated steel
- high quality smooth running plastic housing castors

## warranty

- According to our current service guidelines

medical aid number: (Germany) 22.40.01.0111  
article-ID 8014 manual spreading

# WHAT WERE WE THINKING?

Here are a few USPs of our SLK Universal Hoists!

## Profile



### Dieter Fürst

CAD Designer

A master mechanical engineer and mechanical technician by training, I have been working as a product developer in the field of medical aids / medical products since 2007 and have already developed patient hoists and care beds to market maturity for several manufacturers. Working as a design engineer in automotive manufacturing automation and two years working for a medical device drive manufacturer rounded off my experience until I joined SLK in 2016. This is the right place for development work and innovation. With our team we have already developed several successful hoist models and it is a pleasure to be able to expand our hoist portfolio step by step.

## 1. CENTERING OF THE STANDARD MAST




Patient hoist masts should often be foldable so that the product can be easily transported. At the same time, the mast and chassis should be easily removable without the use of tools. The requirement during use, on the other hand, is for the two parts to be connected as rigidly as possible. To achieve this, we came up with a locking device for the standard mast with a centring piece. This means that the mast is fixed to the chassis without any play. We have protected this invention by means of a German utility model.

## 2. MAXIMUM LOAD 185 KG / 407 LBS



The fact that people are becoming heavier and heavier is also reflected in everyday care. That is why we have developed the SLK Carry hoist family in such a way that we can allow the high maximum load of 185 kg. This covers more than 90% of the weights of patients. Therefore, there is no need to purchase an additional hoist, except for real heavy-duty care.

## 3. TELESCOPIC STANDARD MAST



The telescopic standard mast turns the Multy into a hybrid hoist. It allows the Multy to be used as an active hoist and universal hoist without any restrictions and with an optimal height of the pivot point depending on the type of application. The height difference between the top and bottom position is 340 mm. For an upright raising procedure, the lowest possible pivot point is needed to be able to raise persons with residual mobility upright, while for a transfer using a universal hoist, the highest possible pivot point is needed.

# These are the reasons for using Hybrid Hoists

The SLK Multy is a modular hybrid hoist that can be set up in two ways: As an Active Hoist and as a Universal Hoist. This is done by attaching an appropriate lifting arm and optionally using a foot-board with a shin support. The idea behind this is to provide our customers with a single hoist that can be used as a very well-functioning active hoist while the patient is still mobile. As required or if the patient's mobility decreases later on, the same hoist can be converted into a fully-fledged Universal Hoist / Patient Hoist in just a few steps. The SLK Multy Conversion Kit Active and the SLK Multy Conversion Kit Universal are available for this purpose. This eliminates the need to purchase a second hoist or to purchase a new hoist.

## The advantages of a Hybrid Hoist

### The advantages for the patient

- Quick conversion from an active hoist to a universal hoist and vice versa
- Safe and careful transfer
- Active Hoists provide support in therapy, in maintaining or increasing mobility and independence
- With their features, active hoists are beneficial for the stabilisation of the cardiovascular system
- Active Hoists provide an improvement in body sensation and bowel function
- Active Hoists positively influence the strengthening of bones, joints and cartilage and the careful lifting and setting down of the patient
- Picking up people from the ground possible, electrically (optional) adjustable chassis spreading

### The advantages for the caregiver

- Simple operation during the transfer
- Minimization of the physical strain on the care assistant during the transfer

hybrid hoists  
**designed for use in facilities,  
clinics and homecare**



# SLK Multy Universal / Multy Universal e

Universal Hoist convertible to Active Hoist



SLK Multy Conversion Kit Active

DATA	MANUAL SPREADING	ELECTRIC SPREADING
length	118 cm (46.4")	118 cm (46.4")
total height	137 cm (54")	137 cm (54")
width	61 cm (24")	61 cm (24")
leg height	11 cm (4.3")	11 cm (4.3")
lifting height	81 - 203 cm (32" - 80")	81 - 203 cm (32" - 80")
turning radius	132 cm (52")	132 cm (52")
base expansion	51 - 94 cm (20" - 37")	51 - 94 cm (20" - 37")
weight	37 kg (81 lbs)	38 kg (84 lbs)

Further details on page: 22-23

## facts

- max. patient weight: 185 kg (407 lbs)
- protects the back of the caregiver
- gentle pick up and settle of the patient
- lifting from floor level possible
- electric spreading (optional) adjustable
- high manoeuvrability in a small area
- telescopic stand mast
- Hook mount with spring safety device
- 360-degree rotation of the Spreader Bar in highest possible position
- comfort-control with display and Readout function
- handling via control box possible
- soft-start / soft-stop
- high battery power
- foldable without tools / transportable
- standard mast and chassis are removable without tools

## appropriate slings

- all SLK Slings (from page 25)

## standard

- SLK 2-Point Spreader Bar

## optional

- SLK 4-Point Spreader Bar X: 8060
- SLK 4-Point Spreader Bar H: 8064
- SLK Multy Conversion Kit Active: 8011

Accessories from page 46

## material

- powder coated steel
- ball-bearing, smooth-running comfort-castors with high-quality plastic housing

## warranty

- According to our current service guidelines

**ALSO AVAILABLE WITH e-CHASSIS...**

medical aid number: (Germany) 22.40.01.0100 (manuell), 22.40.01.0112 (electric)  
article-ID 8009 manual spreading, 8013 electric spreading

# SLK Multy Conversion Kit Active

for SLK Multy Universal / Multy Universal e



2 IN 1

MADE  
IN  
GERMANY

LBS  
407



## SLK Multy Conversion Kit Active

The base of the Hybrid Hoist is the same for both versions.

By using attachments (Conversion Kits) the hoist can be converted to an Universal Hoist or an Active Hoist.

This Conversion Kit transforms the SLK Multy Hybrid Hoist (Universal Hoist) into a Active Hoist.

## SLK Multy Conversion Kit Active

- footboard
- shin pads
- Active lifting arm



Height adjustable mast!  
For setting the optimal  
pivot point of both hoist types!

article-ID 8011




# SLK Multy Active / Multy Active e

Active Hoist convertible to Universal Hoist



MADE  
IN  
GERMANY

LBS  
407



SLK Multy Universal Conversion Kit

DATA	MANUAL SPREADING	ELECTRIC SPREADING
length	118 cm (46.4")	118 cm (46.4")
total height	105 cm (42")	105 cm (42")
width	61 cm (24")	61 cm (24")
leg height	11 cm (4.3")	11 cm (4.3")
lifting height	85 - 163 cm (33" - 64")	85 - 163 cm (33" - 64")
turning radius	132 cm (51.5")	132 cm (51.5")
base expansion	51 - 94 cm (20" - 37")	51 - 94 cm (20" - 37")
weight	40 kg (88 lbs)	41 kg (90 lbs)

further details on page: 22-23

## facts

- max. patient weight: 185 kg (407 lbs)
- protects the back of the caregiver
- easier standing up due to particularly deep footboard
- easy transfer through active belt with seat swing/butt support (optional)
- safe transfer even in standing position
- Therapy-supporting, for the preservation or increase mobility
- electric spreading (optional) adjustable
- high manoeuvrability in a small area
- telescopic stand mast
- height adjustable shin pad
- comfort-control with display and Readout function
- handling via control box possible
- soft-start / soft-stop
- high battery power
- foldable without tools / transportable
- standard mast and chassis are removable without tools

## appropriate slings

- all SLK Active Slings (from page 40)

## standard

- SLK Multy shin pad

## optional

- SLK Narrow Calf Band: 8094
  - SLK Multy Conversion Kit Universal: 8010
- Accessories from page 46

## material

- powder coated steel
- ball-bearing, smooth-running comfort-castors with high-quality plastic housing

## warranty

- According to our current service guidelines

**ALSO AVAILABLE WITH e-CHASSIS...**

medical aid number: (Germany) 22.40.01.1001 (manual), 22.40.01.1027 (electric)  
article-ID 8008 manual spreading, 8012 electric spreading



# SLK Multy Conversion Kit Universal

for SLK Multy Active / Multy Active e



2 IN 1

MADE  
IN  
GERMANY

LBS  
407



## SLK Multy Conversion Kit Universal

The base of the Hybrid Hoist is the same for both versions.

By using attachments (Conversion Kits) the hoist can be converted to an Universal Hoist or an Active Hoist.

This Conversion Kit transforms the SLK Multy Hybrid Hoist (Active Hoist) into a Universal Hoist.

## SLK Multy Conversion Kit Universal

- Universal lifting arm
- SLK 2-Point Spreader Bar



Height adjustable mast!  
For setting the optimal  
pivot point of both hoist types!

article-ID 8010

# These are the reasons for using Active Hoists

Active Hoists make a valuable contribution to being able to change positions quickly and easily in the event of incipient mobility impairment. This also helps to preserve the dignity of the patient when discreetly using the toilet. The prerequisite for the use of active hoists is sufficient residual mobility of the patient and the use of the correct active sling.

## The advantages for the patient

- Safe and careful transfer
- Active Hoists provide support in therapy, in maintaining or increasing mobility and independence
- With their features, active hoists are beneficial for the stabilisation of the cardiovascular system
- Active Hoists provide an improvement in body sensation and bowel function
- Active Hoists positively influence the strengthening of bones, joints and cartilage

## The advantages for the caregiver

- Simple operation during the transfer
- Minimization of the physical strain on the caregiver during the transfer

active hoists  
**designed for use in facilities,  
clinics and homecare**



# SLK Eazy-up

The Active Hoist with fixed, particularly low chassis for driving under the seat



ball-bearing dual-casters



optional Calf Band

DATA	FIXED CHASSIS	MANUAL SPREADING
length	99 cm (39")	112 cm (44")
total height	114 cm (45")	115 cm (45.2")
width	58 cm (22")	58 cm (22")
leg height	6,6 cm (2.5")	9 cm (3.5")
lifting height	90 - 159 cm (35.4" - 62.5")	90 - 159 cm (35.4" - 62.5")
turning radius	90 cm (35.4")	114 cm (45")
base expansion		46 - 78 cm (18" - 31")
weight	34 kg (75 lbs)	37 kg (81 lbs)

## facts

- max. patient weight: 160 kg (352 lbs)
- easy to use in narrow spaces and doors
- high manoeuvrability in a small area
- protects the back of the caregiver
- easier standing up due to particularly deep footboard
- Adapted to the natural get up movement
- safe transfer even in standing position
- Therapy-supporting, for the preservation or increase mobility
- pluggable mast
- standard mast and chassis are removable without tools

## appropriate slings

- SLK Active Sling Eco (size M/L)
- SLK Active Sling (size S - XXL)
- SLK Special Active Sling (size S - L)
- SLK Special Active Sling with Seat Swing (size S - XXL)

## standard

- SLK Shin Pad

## optional

- SLK Narrow Calf Band: 8094

## material

- powder coated steel
- ball-bearing dual-casters

## warranty

- According to our current service guidelines


**medical aid number: (Germany) 22.40.01.1050 (fixed chassis), 22.40.01.1049 (manual spreading)**  
**article-ID 8017 fixed chassis, 8016 manual spreading**

technical specifications  
**universal hoists / active hoists**





# Technical Specifications

## Universal Hoists / Active Hoists



	length a	width b	total height c
<b>SLK Carry 185 Eco</b>			
manual spreading	118 (46.4")	61 (24")	137 (54")
<b>SLK Carry 185 Classic / e</b>			
manual spreading	118 (46.4")	61 (24")	137 (54")
electric spreading	118 (46.4")	61 (24")	137 (54")
<b>SLK Carry 185 Pro / e</b>			
manual spreading	115 (45.3")	64 (26")	137 (54")
electric spreading	115 (45.3")	64 (26")	137 (54")
<b>SLK Carry 185 Pro L / e</b>			
manual spreading	115 (45.3")	64 (26")	135 (53")
electric spreading	115 (45.3")	64 (26")	135 (53")
<b>SLK Multy Universal / e</b>			
manual spreading	118 (46.4")	61 (24")	137 (54")
electric spreading	118 (46.4")	61 (24")	137 (54")
<b>SLK Carry Compact</b>			
manual spreading	103 (40.5")	56 (22")	124 (50")
<b>SLK Multy Active / e</b>			
manual spreading	118 (46.4")	61 (24")	105 (42")
electric spreading	118 (46.4")	61 (24")	105 (42")
<b>SLK Eazy-up fix / flex</b>			
fixed chassis	99 (39")	58 (22")	114 (45")
manual spreading	112 (44")	58 (22")	115 (45.2")

\*image DIN EN ISO 10535:2007-04



lifting height on top <b>k*</b>	height of feet <b>d</b>	lifting height below <b>i*</b>	base expansion <b>f</b>	turning radius <b>g</b>	total weight <b>h</b>	max. patient weight <b>i</b>
203 (80")	11 (4.3")	81 (32")	51 - 94 (20"-37")	132 (52")	37 kg (81 lbs)	185 kg (407 lbs)
203 (80")	11 (4.3")	81 (32")	51 - 94 (20"-37")	132 (52")	37 kg (81 lbs)	185 kg (407 lbs)
203 (80")	11 (4.3")	81 (32")	51 - 94 (20"-37")	132 (52")	38 kg (84 lbs)	185 kg (407 lbs)
203 (80")	11 (4.3")	81 (32")	51 - 94 (20"-37")	130 (51")	38 kg (84 lbs)	185 kg (407 lbs)
203 (80")	11 (4.3")	81 (32")	51 - 94 (20"-37")	130 (51")	39 kg (86 lbs)	185 kg (407 lbs)
201 (79")	9 (3.5")	79 (31")	51 - 94 (20"-37")	130 (51")	38 kg (84 lbs)	185 kg (407 lbs)
201 (79")	9 (3.5")	79 (31")	51 - 94 (20"-37")	130 (51")	39 kg (86 lbs)	185 kg (407 lbs)
203 (80")	11 (4.3")	81 (32")	51 - 94 (20"-37")	132 (52")	37 kg (81 lbs)	185 kg (407 lbs)
203 (80")	11 (4.3")	81 (32")	51 - 94 (20"-37")	132 (52")	38 kg (84 lbs)	185 kg (407 lbs)
175 (69")	8,5 (3.4")	81 (32")	45 - 78 (17.7"-30.7")	115 (45")	33 kg (73 lbs)	160 kg (352 lbs)
163 (64")	11 (4.3")	85 (33")	51 - 94 (20"-37")	132 (52")	40 kg (88 lbs)	185 kg (407 lbs)
163 (64")	11 (4.3")	85 (33")	51 - 94 (20"-37")	132 (52")	41 kg (90 lbs)	185 kg (407 lbs)
159 (35.4")	6,6 (2.5")	90 (35.4")		90 (35.4")	34 kg (75 lbs)	160 kg (352 lbs)
159 (35.4")	9 (3.5")	90 (35.4")	46 - 78 (18"-31")	114 (45")	37 kg (81 lbs)	160 kg (352 lbs)

# Quality, Comfort and Safety are our priorities

In the development of our products, special emphasis is placed on therapeutic effectiveness and user-friendliness. For this reason, our hoist slings are specially tailored to the individual needs of patients.

We offer individual solution concepts for children, adults, light or heavy patients, in the homecare sector, in nursing institutions or clinics.

The SLK hoist slings are characterised by their first-class quality and user-friendliness. High-quality materials and professional workmanship ensure a high load capacity of up to 320 kg. In addition, the belts are washable up to a temperature of 70 °C and thus comply with the generally applicable hygiene standards.


The sling loops are colour-coded. This colour coding makes it possible to adjust the hoist sling to the patient more quickly and safely. The different coloured sling edges ensure quick recognition of the sizes. All SLK hoist slings are also suitable for use with hoists from other manufacturers in accordance with the current compatibility declaration. You can find the download for this on our website.

**Download the current compatibility list from our website under Download.**

**Download! [www.slk-gmbh.de](http://www.slk-gmbh.de)**



slings  
**universal slings / active slings**



# SLK Universal Sling

Universal Sling for general application



## facts

- max. patient weight: 320 kg (705 lbs)
- for safe transfer of the patient and lifting from ground level
- ensures a pleasant skin climate during continuous use
- also available in children's sizes (from size XXS)

## sizes / article ID

KIDS		KIDS	
XXS	XS	S	
5221	5222	5113	
M	L	XL	
5114	5115	5116	
XXL	3XL	4XL	
5117	5208	5209	

## Universal Hoists

- all SLK Universal Hoists
- also suitable for hoists of other producers according to our current declaration of conformity (download from our website)

## optional

- SLK Velcro Headrest: 5118
- SLK Extension Loops: 5719

## material / service

- Polyester-Jersey
- washable up to 70 °C (158 °F)


## warranty

- According to our current service guidelines

settlement item no: 22.40.04.0001

# SLK Universal Sling with Headrest

Universal Sling for general application with Headrest



## facts

- max. patient weight: 320 kg (705 lbs)
- snug Universal Sling with integrated Headrest to stabilize the head
- ensures a pleasant skin climate during continuous use
- also available in children's sizes (from size XXS)

## sizes / article ID

KIDS		KIDS
XXS	XS	S
5225	5226	5119
M	L	XL
5120	5121	5122
XXL		
5123		

## Universal Hoists

- all SLK Universal Hoists
- also suitable for hoists of other producers according to our current declaration of conformity (download from our website)

## optional

- SLK Extension Loops: 5719

## material / service

- Polyester-Jersey
- washable up to 70 °C (158 °F)


## warranty

- According to our current service guidelines

settlement item no: 22.40.04.0001

# SLK Universal Bathroom Sling

## Bathing Sling



### facts

- max. patient weight: 320 kg (705 lbs)
- fast drying material
- a net made of polyester jersey enables fast and easy drainage of water
- ensures a pleasant skin climate during continuous use
- also available in children's sizes (from size XXS)

### sizes / article ID

KIDS		KIDS		
XXS	XS	S		
5229	5230	5124		
M	L	XL		
5125	5126	5127		
XXL	3XL	4XL		
5128	5210	5211		

### Universal Hoists

- all SLK Universal Hoists
- also suitable for hoists of other producers according to our current declaration of conformity (download from our website)

### optional

- SLK Velcro Headrest: 5118
- SLK Extension Loops: 5719

### material / service

- Polyester-Jersey
- washable up to 70 °C (158 °F)

### warranty


- According to our current service guidelines

settlement item no: 22.40.04.0001



# SLK Universal Bathroom Sling with Headrest

Universal Sling for bathing with Headrest



## facts

- max. patient weight: 320 kg (705 lbs)
- fast drying material
- the polyester-jersey-net allows a fast and easy water run-off
- ensures a pleasant skin climate during continuous use
- sling with integrated headrest for head stabilization
- also available in children's sizes (from size XXS)

## sizes / article ID

KIDS		KIDS
XXS	XS	S
5227	5228	5129
M	L	XL
5130	5131	5132
XXL		
5133		

## Universal Hoists

- all SLK Universal Hoists
- also suitable for hoists of other producers according to our current declaration of conformity (download from our website)

## optional

- SLK Extension Loops: 5719

## material / service

- Polyester-Jersey
- washable up to 70 °C (158 °F)

## warranty

- According to our current service guidelines

settlement item no: 22.40.04.0001

# SLK Quick Transfer and Toilet Sling

Transfer and Toilet Sling for quick fastening and repositioning



## facts

- max. patient weight: 320 kg (705 lbs)
- easy and quick to fasten and loose
- with back support and chest loop
- with larger cutout for the buttocks
- ensures a pleasant skin climate during continuous use
- also available in children's sizes (from size XXS)

## sizes / article ID

KIDS		KIDS		
XXS	XS	S		
5219	5220	5134		
M	L	XL		
5135	5136	5137		

## Universal Hoists

- all SLK Universal Hoists
- also suitable for hoists of other producers according to our current declaration of conformity (download from our website)

## optional

- SLK Adjustable Headrest Reinforced: 5138
- SLK Extension Loops: 5719

## material / service

- Polyester-Jersey
- washable up to 70 °C (158 °F)


## warranty

- According to our current service guidelines

settlement item no: 22.40.04.0001

# SLK Toilet Sling

## Toilet Sling



**LBS**  
705

### facts

- max. patient weight: 320 kg (705 lbs)
- quick to put on and put off
- with larger cutout for the buttocks
- ensures a pleasant skin climate during continuous use

### sizes / article ID

<b>S</b> 5139	<b>M</b> 5140	<b>L</b> 5141
<b>XL</b> 5142	<b>XXL</b> 5143	<b>3XL</b> 5206
<b>4XL</b> 5207		

### Universal Hoists

- all SLK Universal Hoists
- also suitable for hoists of other producers according to our current declaration of conformity (download from our website)

### optional

- SLK Adjustable Headrest Reinforced: 5138
- SLK Extension Loops: 5719

### material / service

- Polyester-Jersey
- washable up to 70 °C (158 °F)

### warranty


- According to our current service guidelines

settlement item no: 22.40.04.0001



# SLK Hammock Universal Sling

Hammock Cut Sling for general application



## facts

- max. patient weight: 320 kg (705 lbs)
- high pressure relief for the patient
- specially developed for leg amputees
- ensures a pleasant skin climate during continuous use

## sizes / article ID

<b>S</b> 5144	<b>M</b> 5145	<b>L</b> 5146
<b>XL</b> 5147	<b>XXL</b> 5148	<b>3XL</b> 5200
<b>4XL</b> 5201		

## Universal Hoists

- all SLK Universal Hoists
- also suitable for hoists of other producers according to our current declaration of conformity (download from our website)

## optional

- SLK Velcro Headrest: 5118
- SLK Extension Loops: 5719

## material / service

- Polyester-Jersey
- washable up to 70 °C (158 °F)

## warranty


- According to our current service guidelines

settlement item no: 22.40.04.0001



# SLK Hammock Bathroom Sling

Hammock Cut Sling for general application in the bathroom



## facts

- max. patient weight: 320 kg (705 lbs)
- high pressure relief for the patient
- specially developed for leg amputees
- for transfer to bathtub, shower and pool
- ensures a pleasant skin climate during continuous use

## sizes / article ID

<b>S</b> 5149	<b>M</b> 5150	<b>L</b> 5151
<b>XL</b> 5152	<b>XXL</b> 5153	<b>3XL</b> 5202
<b>4XL</b> 5203		

## Universal Hoists

- all SLK Universal Hoists
- also suitable for hoists of other producers according to our current declaration of conformity (download from our website)

## optional

- SLK Velcro Headrest: 5118
- SLK Extension Loops: 5719

## material / service

- Polyester-Jersey
- washable up to 70 °C (158 °F)


## warranty

- According to our current service guidelines

settlement item no: 22.40.04.0001

# SLK Hammock Seat Pan Sling

Seat Pan Sling for remaining in the wheelchair



**LBS**  
705

## facts

- max. patient weight: 320 kg (705 lbs)
- supple fabric, adjusting to the contours of the seat shell
- no seams or loops in the back and buttocks area
- ensures a pleasant skin climate during continuous use
- also available in children's sizes (from size XXS)

## sizes / article ID

KIDS		
XXS	XS	S
5223	5224	5154
M	L	XL
5155	5156	5157
XXL	3XL	4XL
5158	5204	5205

## Universal Hoists

- all SLK Universal Hoists
- also suitable for hoists of other producers according to our current declaration of conformity (download from our website)

## optional

- SLK Velcro Headrest: 5118
- SLK Extension Loops: 5719

## material / service

- Polyester-Jersey
- washable up to 70 °C (158 °F)


## warranty

- According to our current service guidelines

settlement item no: 22.40.04.0001

# SLK Comfort Lying Sling

Lying Transfer Sling for almost lying transport



## facts

- max. patient weight: 320 kg (705 lbs)
- application with adjustable headrest recommended
- ensures a pleasant skin climate during continuous use

## sizes / article ID

<b>S</b> 5159	<b>M</b> 5160	<b>L</b> 5161
<b>XL</b> 5162	<b>XXL</b> 5163	

## Universal Hoists

- all SLK Universal Hoists
- also suitable for hoists of other producers according to our current declaration of conformity (download from our website)

## optional

- SLK Velcro Headrest: 5118
- SLK Extension Loops: 5719

## material / service

- Polyester-Jersey
- washable up to 70 °C (158 °F)

## warranty

- According to our current service guidelines

settlement item no: 22.40.04.0001



# SLK Comfort Lying Bathroom Sling

Lying Transfer Sling for almost lying transport



## facts

- max. patient weight: 320 kg (705 lbs)
- high pressure relief for the patient
- specially developed for leg amputees
- for transfer to bathtub, shower and pool
- ensures a pleasant skin climate during continuous use

## sizes / article ID

<b>S</b> 5164	<b>M</b> 5165	<b>L</b> 5166
<b>XL</b> 5167	<b>XXL</b> 5168	

## Universal Hoists

- all SLK Universal Hoists
- also suitable for hoists of other producers according to our current declaration of conformity (download from our website)

## optional

- SLK Velcro Headrest: 5118
- SLK Extension Loops: 5719

## material / service

- Polyester-Jersey
- washable up to 70 °C (158 °F)


## warranty

- According to our current service guidelines

settlement item no: 22.40.04.0001

# SLK Back-Fix Sling

## Universal Sling



### facts

- max. patient weight: 320 kg (705 lbs)
- with reinforced back and integrated headrest
- ensures a pleasant skin climate during continuous use
- for patients whose back and head need to be extremely supported

### sizes / article ID

S  
5169

M  
5170

L  
5171

XL  
5172

### Universal Hoists

- all SLK Universal Hoists
- also suitable for hoists of other producers according to our current declaration of conformity (download from our website)

### optional

- SLK Extension Loops: 5719

### material / service

- Polyester-Jersey
- washable up to 70 °C (158 °F)

### warranty

- According to our current service guidelines

settlement item no: 22.40.04.0001

# SLK Back-Fix Bathroom Sling

## Bathing Sling



### facts

- max. patient weight: 320 kg (705 lbs)
- with reinforced back and integrated headrest
- ensures a pleasant skin climate during continuous use
- for transfer to bathtub, shower and pool
- for patients whose back and head need to be extremely supported

### sizes / article ID

S 5173	M 5174	L 5175
XL 5176		

### Universal Hoists

- all SLK Universal Hoists
- also suitable for hoists of other producers according to our current declaration of conformity (download from our website)

### optional

- SLK Extension Loops: 5719

### material / service

- Polyester-Jersey
- washable up to 70 °C (158 °F)

### warranty

- According to our current service guidelines

settlement item no: 22.40.04.0001



# SLK Maxi Sling

Universal Sling for bariatric patients



## facts

- max. patient weight: 320 kg (705 lbs)
- with integrated headrest
- ensures a pleasant skin climate during continuous use
- extra large sizes

## sizes / article ID

XXL  
5177

3XL  
5178

4XL  
4XL

## Universal Hoists

- all SLK Universal Hoists
- also suitable for hoists of other producers according to our current declaration of conformity (download from our website)

## optional

- SLK Extension Loops: 5719

## material / service

- Polyester-Jersey
- washable up to 70 °C (158 °F)

## warranty

- According to our current service guidelines

settlement item no: 22.40.04.0001

# SLK Active Sling

Raising Sling for dual- and upward pulling Active Hoists



## facts

- max. patient weight: 320 kg (705 lbs)
- padded sling supporting raising motion
- additional abdominal belt for greater safety
- high wearing comfort

## sizes / article ID

<b>S</b> 5100	<b>M</b> 5101	<b>L</b> 5102
<b>XL</b> 5103	<b>XXL</b> 5104	

## Active Hoists

- all SLK Active Hoists
- also suitable for hoists of other producers according to our current declaration of conformity (download from our website)

## optional

- SLK Extension Loops: 5719

## material / service

- Polyester-Jersey
- washable up to 70 °C (158 °F)

## warranty

- According to our current service guidelines

settlement item no: 22.40.04.0001



# SLK Special Active Sling

Raising Sling with buttocks support for upward pulling Active Hoists



## facts

- max. patient weight: 320 kg (705 lbs)
- padded sling ensuring a safe and comfortable feeling while getting up
- additional buttocks belt for greater safety

## sizes / article ID

S	M	L
5105	5106	5107

## Active Hoists

- all SLK Active Hoists
- also suitable for hoists of other producers according to our current declaration of conformity (download from our website)

## optional

- SLK Extension Loops: 5719

## material / service

- Polyester-Jersey
- washable up to 70 °C (158 °F)


## warranty

- According to our current service guidelines

settlement item no: 22.40.04.0001

# SLK Special Active Sling with Seat Swing

Raising Sling for upward pulling Active Hoists



**LBS**  
**705**

## facts

- max. patient weight: 320 kg (705 lbs)
- padded sling ensuring a safe and comfortable feeling while getting up

## sizes / article ID

<b>S</b> 5188	<b>M</b> 5189	<b>L</b> 5190
<b>XL</b> 5191	<b>XXL</b> 5192	

## Active Hoists

- all SLK Active Hoists
- also suitable for hoists of other producers according to our current declaration of conformity (download from our website)

## optional

- SLK Extension Loops: 5719

## material / service


- Polyester-Jersey
- washable up to 70 °C (158 °F)

## warranty

- According to our current service guidelines

settlement item no: 22.40.04.0001

slings  
size chart / overview




## Body height in cm


Weight in kg

	80	90	100	110	120	130	140	150	160	170	180	190	200	210
15	KIDS XXS	KIDS XXS	KIDS XXS	KIDS XXS										
20	KIDS XXS	KIDS XXS	KIDS XS	KIDS XS	KIDS XS	KIDS XS								
25	KIDS XS	KIDS XS	KIDS XS	KIDS XS	KIDS XS	KIDS XS	KIDS XS							
30			S	S	S	S	S	S	M	M				
40			S	S	S	S	S	S	M	M	M			
50			M	M	M	M	S	S	M	M	M	M	M	M
60			M	M	M	M	M	M	M	M	M	M	M	M
70			M	M	M	M	M	M	L	L	L	M	M	M
80				L	L	L	L	L	L	L	L	L	L	L
90				L	L	L	L	L	L	L	L	L	L	L
100								L	L	L	L	L	L	L
110								XL	XL	XL	XL	L	L	L
120								XXL	XL	XL	XL	XL	XL	XL
130								XXL	XL	XL	XL	XL	XL	XL
140								XXL	XL	XL	XL	XL	XL	XL
150								XXL	XXL	XXL	XXL	XXL	XXL	XXL
160								XXL	XXL	XXL	XXL	XXL	XXL	XXL
170								XXL	XXL	XXL	XXL	XXL	XXL	XXL
180								XXL	XXL	XXL	XXL	XXL	XXL	XXL
190									XXL	XXL	XXL	XXL	XXL	XXL
200									XXL	XXL	XXL	XXL	XXL	XXL
210									XXL	XXL	XXL	3XL	3XL	3XL
220									3XL	3XL	3XL	3XL	3XL	3XL
230									3XL	3XL	3XL	3XL	3XL	3XL
240									3XL	3XL	3XL	3XL	3XL	3XL
250									3XL	3XL	3XL	3XL	3XL	3XL
260									3XL	4XL	4XL	3XL	3XL	3XL
270									3XL	4XL	4XL	4XL	4XL	4XL
280										4XL	4XL	4XL	4XL	4XL
290										4XL	4XL	4XL	4XL	4XL
300										4XL	4XL	4XL	4XL	4XL
310										4XL	4XL	4XL	4XL	4XL
320										4XL	4XL	4XL	4XL	4XL


SLK Velcro Headrest



SLK Adjustable Headrest  
Reinforced



SLK Toilet Sling



SLK Calf Band




SLK Special Active Sling



SLK Special Active Sling  
with Seat Swing



SLK Quick Transfer and Toilet Sling




SLK Active Sling




SLK Seat Swing




SLK Universal Sling




SLK Universal Sling  
with Headrest




SLK Back-Fix Sling




SLK Maxi Sling




SLK Hammock Seat  
Pan Sling




SLK Hammock  
Universal Sling



SLK Comfort Lying Sling



accessories  
**hoist accessories**



# SLK Hoist Scales

## Hoist Scales for all SLK Universal Hoists



### SLK Hoist Scale

The SLK hoist scale can be connected directly to the universal hoist and enables people to be weighed.

It is easy and quick to integrate into the daily care routine. The patient scale is designed for a carrying weight of up to 300 kg and therefore suitable for bariatric patients.

- including adapter for all SLK Universal Hoists
- Use of hoists other manufacturers on request

**article-ID 5186**



### SLK Calibrated Hoist Scale

The SLK verifiable hoist scale is factory-verified and can therefore be used for medical purposes. It can be connected directly to the patient hoist and enables people to be weighed.

It can be easily and quickly integrated into the daily care routine. The patient scale is designed for a carrying weight of up to 300 kg and is therefore also suitable for bariatric patients.

- including adapter for all SLK Universal Hoists
- Use of hoists other manufacturers on request

**article-ID 5183**

# SLK Carry Charge Base

## SLK Charger Propulsion system Typ I + Typ II



With the SLK Carry charging station batteries of patient hoists can be charged independently of the location of the hoist. This can be used as a wall charging station or as a table charging station. Use in stationary facilities is recommended, as with an additional exchangeable battery the hoist is ready for use at any time.

- suitable for SLK Carry Productfamily, SLK Multy Hybrid Hoist
- Charger optionally available
- Delivery without battery

**Typ I: article-ID 8634**  
**Typ II: article-ID 8454**

## SLK Charging Station with Charger Propulsion system Typ I + Typ II



With the SLK Carry charging station batteries of patient hoists can be charged independently of the location of the hoist. This can be used as a wall charging station or as a table charging station. Use in stationary facilities is recommended, as with an additional exchangeable battery the hoist is ready for use at any time.

- suitable for SLK Carry Productfamily, SLK Multy Hybrid Hoist
- Delivery without battery

**Typ I: article-ID 8641**  
**Typ II: article-ID 8390**

## SLK Wall-mounted Charging Station UMEC Propulsion system



With the SLK Carry charging station batteries of patient hoists can be charged independently of the location of the hoist. This can be used as a wall charging station or as a table charging station. Use in stationary facilities is recommended, as with an additional exchangeable battery the hoist is ready for use at any time.

- suitable for all SLK Universal Hoists (except SLK Carry's + Multy's)
- Delivery without battery

**article-ID 5182**



# SLK Spreader Bar

## SLK 2-Point and 4-Point Spreader Bar



■ width: 59 cm (23.2")



### SLK 2-Point Spreader Bar

The SLK 2-Point Spreader Bar allows a particularly comfortable transport feeling. Slings can not be unintentionally released due to its spring lock.

#### article-ID 8053 (telegrey)

■ for SLK Carry 185 Classic / e

#### article-ID 8054 (silver)

■ for SLK Carry 185 Pro, Pro L, SLK Multy Universal and the electrically spreadable models

#### article-ID 8689 (telegrey)

■ for SLK Carry Compact Universal Hoist



■ width: 60 cm (23.6")



### SLK 4-Point Spreader Bar X

The SLK 4-Point Spreader Bar X allows a particularly comfortable transport feeling. Slings can not be unintentionally released due to its spring lock.

#### article-ID 8060 (silver)

■ for SLK Carry Productfamily, SLK Multy Universal and the electrically spreadable models



■ width: 57 cm (22.4")



### SLK 4-Point Spreader Bar H

The SLK 4-Point Spreader Bar H allows a particularly comfortable transport feeling. It is especially suitable for obese patients. Slings can not be unintentionally released due to its spring lock.

**We recommend the use of the SLK 4-Point Spreader Bar H inter alia to:**

- obese patients
- sensitive patients
- patients which tend to stretch spastik



■ width: 52 cm (20.5")



#### article-ID 8064 (silver)

■ for SLK Carry Productfamily, SLK Multy Universal and the electrically spreadable models

#### article-ID 8095 (telegrey)

■ specially for SLK Carry Compact Universal Hoist

# SLK 8-Point Spreader Bar

Tareable 8-Point Spreader Bar for gentle Horizontal transfer



SLK 8-Point Lying Sling



Hook mount with spring lock



The purpose of the SLK 8-Point Spreader Bar is to enable lifting and transferring of persons together with a lift system and an 8-point sling. The special feature of the SLK 8-Point Spreader Bar is that this transfer of patients is performed in the lying position.

As suitable Universal Hoist systems we name the SLK Carry and Multy Universal Hoist.

## facts

- Tareable 8-Point Spreader Bar for gentle Horizontal transfer
- Ideal for transferring comatose and pain sensitive patients
- Can be used in conjunction with 8-Point Lying Sling

## optional

- SLK Velcro Headrest: 5118

## scope of delivery

- SLK 8-Point Spreader Bar
- 1 x instruction manual 8303

## article-ID

**8062** (Maximum load 185 kg/407 lbs for SLK Carry 185 Productfamily and SLK Multy Universal)

**8063** (Maximum load 275 kg/606 lbs)

**8065** (Maximum load 300 kg/660 lbs for SLK Flyer XL)

**8059** (Maximum load 400 kg/880 lbs)

**5699** (suitable SLK 8-Point Lying Sling)

## technical data / dimensions

- length: 161 cm (63")
- width: 74,5 cm (29")
- height: 29 cm (11")
- weight: 17 kg - 19 kg (37" - 42")  
(depending on version)

# SLK Hoist Accessories



## SLK Impact Protector for Lifting Arm

The SLK Impact Protector for Lifting Arm reduces Injury risks.

- suitable for all SLK Universal Hoists
- can also be mounted on Universal Hoists from other manufacturers

**article-ID 5464**



## SLK Velcro Headrest

The SLK Velcro Headrest is optional for supporting the patient's head and neck area. For ideal SLK 4-Point Spreader Bar is recommended for ideal use.

- suitable for Slings with matching device

**article-ID 5118**



## SLK Adjustable Headrest Reinforced

The Adjustable Headrest Reinforced is optional to support the patient's head, neck and back. For ideal application a SLK 4-Point Spreader Bar is recommended.

- suitable for SLK Toilet Sling and SLK Quick Transfer and Toilet Sling

**article-ID 5138**

# SLK Hoist Accessories



## SLK Extension Loops for the Hoist Slings

The SLK extension loops enable you to extend the belt loops of suitable SLK hoist belts in order to adapt them quickly and easily to the different needs of patients. They are made of extremely tear-resistant webbing and are easy to attach to the webbing loops of SLK lift straps and are intended for SLK lift straps on patient hoists.

**By extending the strap loops, the following goals can be achieved:**

- Increasing the free space of the patient to the transport bar
  - easier hooking of the belt loops on the transport bar
  - especially when picking up the patient from the floor
  - Extension of the leg loops in order to achieve a more upright position or to make the angle of the position more position more open (in the hip)
  - Effective extension, after attachment to SLK Slings approx. 18 cm
- suitable for all SLK Slings
  - standard: DIN EN ISO 10535:2007-04
  - length: 24 cm (folded)

**article-ID 5719**

# #workgreen

# SLK Weight Test Trolley up to 220 kg

for maximal load test as a part of the hoist maintenance



MADE  
IN  
GERMANY

LBS  
485



without load



suspended condition

## facts

- max. load: 220 kg (485 lbs)
- for the functional testing of Universal Hoists with maximum load according to EN 10535 (appendix B1)
- adapt to your patient hoist the height of the stop can be simply determined by the weight trolley, therefore the whole lifting area is covered during the test
- suitable for almost all common hoists

## material / service

- powder coated steel
- braked castors 75 ø mm (3"), ball bearing

## warranty

- According to our current service guidelines

## scope of delivery

- frame with castors **without weights**
- drawbar
- shackles
- bott with sl-fastener

# SLK Weight Test Trolley up to 440 kg

for maximal load test as a part of the hoist maintenance



MADE  
IN  
GERMANY

LBS  
970



suspended condition  
(Test with a 185 kg hoist)



with load

## facts

- max. load: 440 kg (970 lbs)
- for the functional testing of Universal Hoists with maximum load according to EN 10535 (appendix B1)
- adapt to your patient hoist the height of the stop can be simply determined by the weight trolley, therefore the whole lifting area is covered during the test
- suitable for almost all common hoists

## material / service

- powder coated steel
- braked castors 75 ø mm (3"), ball bearing

## warranty

- According to our current service guidelines

## scope of delivery

- 2 frames with castors **without weights**
- drawbar
- 2 shackles
- 2 Slings
- bott with sl-fastener

article-ID 5474





■ MADE  
■ IN  
■ GERMANY





## Products

- alternating pressure systems for decubitus prophylaxis and therapy
- intermittent pneumatic compression with 12 chamber garments for different fields of applications
- hoists for patients with mobility restrictions due to illness or disability

## Service

- fast delivery of products
- safe and timely chemothermic cleaning and testing of anti-decubitus systems according to a validated procedure
- regular inspection of the cleaning by an internationally recognized hygiene institute
- cleaned systems are packed in breathable sanitary bags, sealed and shipped in new boxes
- maintenance and repair of SLK products
- cleanings with maintenance and leakage test according to DIN EN 62353
- repairs of external system on request
- high security for re-use

## Training

- SLK offers advanced product training programs as well as trainings about bedsores for physicians, nursing personnel, hospitals and / or the medical service employees and others
- we offer training for our partners on the topic of wound treatment



## Wellell Group

SLK Vertriebsgesellschaft mbH  
Am Herdicksbach 18  
45731 Waltrop - Germany

telephone: +49 (0) 2309 - 91545 - 0  
telefax: +49 (0) 2309 - 91545 - 999  
e-mail: [info@slk-gmbh.de](mailto:info@slk-gmbh.de)  
internet: [www.slk-gmbh.de](http://www.slk-gmbh.de)  
version: 1.7 - date: july 2022

 @SLK Medizintechnik

 [slk\\_medizintechnik](https://www.instagram.com/slk_medizintechnik)